

GLUT3 Rabbit Monoclonal Antibody



CAB4137

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

48-60kDa

Calculated MW:

54kDa

Applications:

WB

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

6515

Uniprot

P11169

Synonyms:

GLUT3


Immunogen:

A synthesized peptide derived from human GLUT3

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH7.3.

Product Images



Western blot - GLUT3 Rabbit mAb (CAB4137)