

CAB4149

Product Information

Product SKU:	CAB4149	Gene ID:	10411	Size:	20uL, 100uL
Clone No:	ARC0908	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	100kDa	Conjugate:	Unconjugated
Calculated MW:	104kDa	Isotype:	IgG

Immunogen Information

Background: Enables guanyl-nucleotide exchange factor activity and protein domain specific binding activity. Involved in several processes, including positive regulation of protein modification process

Recommended Dilution: WB,1:500 - 1:1000 IHC-P,1:50 - 1:200

Synonyms: EPAC; Epac1/RAPGEF3

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 824-923 of human Epac1/RAPGEF3 (O95398).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.