

PGAM1 Rabbit Monoclonal Antibody



CAB4170

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

28KDa

Calculated MW:

29kDa

Applications:

WB IHC

Reactivity:

Human, Mouse, Rat

Protein Background

The protein encoded by this gene is a mutase that catalyzes the reversible reaction of 3-phosphoglycerate (3-PGA) to 2-phosphoglycerate (2-PGA) in the glycolytic pathway. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Nov 2015]

Immunogen information

Gene ID:

5223

Uniprot

P18669

Synonyms:

HEL-S-35; PGAM-B; PGAMA

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IHC 1:50
- 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification


Immunogen:

A synthesized peptide derived from human PGAM1

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH7.3.

Product Images



Western blot - PGAM1 Rabbit mAb (CAB4170)