

UGT2B4 Rabbit Polyclonal Antibody



CAB4180

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

50kDa

Calculated MW:

42kDa/44kDa/60kDa

Applications:

WB

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:200 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

7363

Uniprot

P06133

Synonyms:

UGT2B4; HLUG25; UDPGT2B4; UDPGTH1; UDPGTh-1; UGT2B11


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 280-480 of human UGT2B4 (NP_066962.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot analysis of extracts of various cell lines, using UGT2B4 antibody (CAB4180) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (CABS014) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (CABM00020). Exposure time: 10s.