

TRPM8 Rabbit Monoclonal Antibody



CAB4269

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

110KDa

Calculated MW:

128kDa

Applications:

WB

Reactivity:

Human

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

79054

Uniprot

Q7Z2W7

Synonyms:

LTRPC6; TRPP8


Immunogen:

A synthesized peptide derived from human TRPM8

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH7.3.

Product Images



Western blot - TRPM8 Rabbit mAb (CAB4269)