

**CAB4286**

## Product Information

<b>Product SKU:</b>	CAB4286	<b>Gene ID:</b>	9588	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	ARC0954	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

## Additional Information

<b>Observed MW:</b>	25kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	25kDa	<b>Isotype:</b>	IgG

## Immunogen Information

**Background:** The protein encoded by this gene is a member of the thiol-specific antioxidant protein family. This protein is a bifunctional enzyme with two distinct active sites. It is involved in redox regulation of the cell; it can reduce H(2)O(2) and short chain organic, fatty acid, and phospholipid hydroperoxides. It may play a role in the regulation of phospholipid turnover as well as in protection against oxidative injury.

**Recommended Dilution:** WB,1:500 - 1:1000

**Synonyms:** PRX; p29; AOP2; 1-Cys; NSGPx; aiPLA2; LPCAT-5; HEL-S-128m; Peroxiredoxin 6

**Purification Method:** Affinity purification

**Immunogen:** A synthetic peptide corresponding to a sequence within amino acids 125-224 of human Peroxiredoxin 6 (NP\_004896.1).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.