

Product Information

Product SKU:	CAB4294	Gene ID:	9373	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	87kDa	Conjugate:	Unconjugated
Calculated MW:	87kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable ubiquitin binding activity. Involved in cellular response to lipopolysaccharide; macroautophagy; and positive regulation of phospholipase A2 activity. Located in cytoplasm; extracellular exosome; and nucleus.

Recommended Dilution:	WB,1:200 - 1:2000
Synonyms:	DOA1; PLAP; PLA2P; NDMSBA; PLAA
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 536-795 of human PLAA (NP_001026859.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.