

## CAB4442

---

### Product Information

<b>Product SKU:</b>	CAB4442	<b>Gene ID:</b>	10625	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

---

### Additional Information

<b>Observed MW:</b>	72kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	72kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** Involved in RNA splicing; negative regulation of protein ubiquitination; and response to virus. Located in cytosol. Implicated in immunodeficiency 70.

**Recommended Dilution:** WB,1:500 - 1:2000

**Synonyms:** ND1; ARA3; NS-1; IMD70; NS1BP; FLARA3; KLHL39; NS1-BP; HSPC068; IVNS1ABP

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human IVNS1ABP (NP\_006460.2).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.