

CAB4894

Product Information

| | | | | | |
|--------------|---------|---------------|--------|-------------|-----------------|
| Product SKU: | CAB4894 | Gene ID: | 9367 | Size: | 20uL, 100uL |
| Clone No: | ARC0332 | Host Species: | Rabbit | Reactivity: | Human,Mouse,Rat |

Additional Information

| | | | |
|----------------|-------|------------|--------------|
| Observed MW: | 23kDa | Conjugate: | Unconjugated |
| Calculated MW: | 23kDa | Isotype: | IgG |

Immunogen Information

Background: Enables GDP binding activity; GTP binding activity; and GTPase activity. Involved in negative regulation by host of symbiont catalytic activity; positive regulation of exocytosis; and regulation of protein localization. Located in late endosome; lysosome; and phagocytic vesicle.

Recommended Dilution: WB,1:500 - 1:2000 IHC-P,1:50 - 1:200 IF/ICC,1:50 - 1:200

Synonyms: RAB9; RAB9A

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human RAB9A (NP_004242.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.