

Product Information

Product SKU:	CAB5064	Gene ID:	8451	Size:	20uL, 100uL
Clone No:	ARC1237	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	88kDa	Conjugate:	Unconjugated
Calculated MW:	88kDa	Isotype:	IgG

Immunogen Information

Background: CUL4A is the ubiquitin ligase component of a multimeric complex involved in the degradation of DNA damage-response proteins (Liu et al., 2009 [PubMed 19481525]).

Recommended Dilution:	WB,1:500 - 1:2000 IHC-P,1:50 - 1:200
Synonyms:	CUL4A; cullin-4A; Cullin 4A
Purification Method:	Affinity purification
Immunogen:	A synthetic peptide corresponding to a sequence within amino acids 660-759 of human Cullin 4A (Q13619).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.