

CAB5069

Product Information

Product SKU:	CAB5069	Gene ID:	84869	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Rat

Additional Information

Observed MW:	25kDa	Conjugate:	Unconjugated
Calculated MW:	25kDa	Isotype:	IgG

Immunogen Information

Background: Enables several functions, including 3-oxoacyl-[acyl-carrier-protein] reductase (NADPH) activity; NADPH binding activity; and NADPH dehydrogenase (quinone) activity. Involved in fatty acid biosynthetic process; glycoside metabolic process; and protein tetramerization. Located in mitochondrial matrix. Part of oxidoreductase complex.

Recommended Dilution: WB,1:200 - 1:2000

Synonyms: SDR45C1; CBR4

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-237 of human CBR4 (NP_116172.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.