

CAB5108

Product Information

Product SKU:	CAB5108	Gene ID:	113251	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse,Rat

Additional Information

Observed MW:	81kDa	Conjugate:	Unconjugated
Calculated MW:	81kDa	Isotype:	IgG

Immunogen Information

Background: Enables mRNA 3'-UTR binding activity and poly(A) binding activity. Involved in cytoskeleton organization; positive regulation of translation; and regulation of cell morphogenesis. Located in cytosol. Colocalizes with cytoplasmic stress granule; cytosolic small ribosomal subunit; and polysome.

Recommended Dilution:	WB,1:500 - 1:1000 IF/ICC,1:50 - 1:200
Synonyms:	PP13296; LARP4
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of human LARP4 (NP_443111.4).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.