

CAB5128

Product Information

Product SKU:	CAB5128	Gene ID:	121441	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	78kDa	Conjugate:	Unconjugated
Calculated MW:	72kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to be involved in protein localization to centrosome. Located in centrosome; nucleoplasm; and plasma membrane.

Recommended Dilution: WB,1:500 - 1:2000 IF/ICC,1:50 - 1:200

Synonyms: GCP-WD; TUBGCP7; NEDD1

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 411-660 of human NEDD1 (NP_690869.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.