

SUN5 Rabbit Polyclonal Antibody



CAB5158

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

39kDa/42kDa

Applications:

WB

Reactivity:

Human, Mouse

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

The protein encoded by this gene appears to play a role in the meiotic stage of spermatogenesis. The encoded protein localizes to the junction between the sperm head and body and may be involved in nuclear envelope reconstitution and nuclear migration. Mutations in this gene have been implicated in acephalic spermatozoa syndrome.

Immunogen information

Gene ID:

76407

Uniprot

Synonyms:

SUN5; SPAG4L; SPGF16; TSARG4; dJ726C3.1

Immunogen:

Recombinant protein of mouse SUN5

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
