

CAB5236

Product Information

Product SKU:	CAB5236	Gene ID:	1999	Size:	20uL, 100uL
Clone No:	ARC1191	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	42kDa	Conjugate:	Unconjugated
Calculated MW:	41kDa	Isotype:	IgG

Immunogen Information

Background: Enables DNA-binding transcription activator activity, RNA polymerase II-specific and sequence-specific double-stranded DNA binding activity. Involved in inflammatory response; negative regulation of transcription, DNA-templated; and positive regulation of transcription by RNA polymerase II. Located in Golgi apparatus; cytosol; and nucleoplasm.

Recommended Dilution: WB,1:500 - 1:2000 ChIP,5µg antibody for 10µg-15µg of ChromatinCUT&Tag, 10⁵ cells /1 µg

Synonyms: ERT; ESX; EPR-1; ESE-1; ESE1/ELF3

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-371 of human ESE1/ELF3 (P78545).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.