

# EIF4A1 Rabbit Polyclonal Antibody



CAB5294

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

45KDa

### Calculated MW:

39kDa/46kDa

### Applications:

WB IF IP

### Reactivity:

Human, Mouse, Rat

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000 IF 1:20 -  
1:50 IP 1:20 - 1:50

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

1973

#### Uniprot

P60842

#### Synonyms:

EIF4A1; DDX2A; EIF-4A; EIF4A; eIF-4A-I; eIF4A-I

#### Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-406 of human EIF4A1 (NP\_001407.1).

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

