

Product Information

Product SKU:	CAB5322	Gene ID:	51588	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	57-75kDa	Conjugate:	Unconjugated
Calculated MW:	57kDa	Isotype:	IgG

Immunogen Information

Background: Enables SUMO ligase activity and ubiquitin protein ligase binding activity. Involved in negative regulation of transcription, DNA-templated; positive regulation of protein sumoylation; and protein sumoylation. Located in cytoplasm and nucleus.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	PIASY; Piasg; ZMIZ6; PIAS-gamma; PIAS4
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 60-260 of human PIAS4 (NP_056981.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.