

Product Information

Product SKU:	CAB5369	Gene ID:	8065	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	91kDa	Conjugate:	Unconjugated
Calculated MW:	91kDa	Isotype:	IgG

Immunogen Information

Background: Enables ubiquitin protein ligase binding activity. Predicted to be involved in SCF-dependent proteasomal ubiquitin-dependent protein catabolic process and protein ubiquitination. Predicted to act upstream of or within cerebral cortex radially oriented cell migration and radial glia guided migration of Purkinje cell. Located in site of DNA damage. Part of Cul5-RING ubiquitin ligase complex.

Recommended Dilution:	WB,1:1000 - 1:5000
Synonyms:	CUL-5; VACM1; VACM-1; CUL5
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human CUL5 (NP_003469.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.