

Product Information

Product SKU:	CAB5381	Gene ID:	5905	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	64kDa	Conjugate:	Unconjugated
Calculated MW:	64kDa	Isotype:	IgG

Immunogen Information

Background:	This gene encodes a protein that associates with the nuclear pore complex and participates in the regulation of nuclear transport. The encoded protein interacts with Ras-related nuclear protein 1 (RAN) and regulates guanosine triphosphate (GTP)-binding and exchange. Alternative splicing results in multiple transcript variants.
Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	SD; Fug1; RANGAP; RanGAP1
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 398-587 of human RanGAP1 (NP_002874.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.