

---

## Product Information

<b>Product SKU:</b>	CAB5501	<b>Gene ID:</b>	7088	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human

---

## Additional Information

<b>Observed MW:</b>	83kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	83kDa	<b>Isotype:</b>	IgG

---

## Immunogen Information

**Background:** Enables identical protein binding activity and transcription corepressor activity. Involved in negative regulation of I-kappaB kinase/NF-kappaB signaling; negative regulation of anoikis; and regulation of gene expression. Located in cytosol and nucleoplasm. Part of beta-catenin-TCF complex.

**Recommended Dilution:** WB, 1:500 - 1:2000

**Synonyms:** ESG; ESG1; GRG1; TLE-1; TLE1

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 501-770 of human TLE1 (NP\_005068.2).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.