

Product Information

Product SKU:	CAB5550	Gene ID:	51343	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	55kDa	Conjugate:	Unconjugated
Calculated MW:	55kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable anaphase-promoting complex binding activity and ubiquitin ligase activator activity. Involved in anaphase-promoting complex-dependent catabolic process; mitotic G2 DNA damage checkpoint signaling; and positive regulation of protein metabolic process. Located in nuclear membrane and nucleoplasm. Colocalizes with anaphase-promoting complex.

Recommended Dilution:	WB,1:100 - 1:500 IF/ICC,1:50 - 1:200
Synonyms:	FZR; CDH1; FZR2; HCDH; HCDH1; CDC20C; DEE109; FZR1
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 224-493 of human FZR1 (NP_057347.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.