

CAB5599

Product Information

Product SKU:	CAB5599	Gene ID:	117581	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse,Rat

Additional Information

Observed MW:	18kDa	Conjugate:	Unconjugated
Calculated MW:	18kDa	Isotype:	IgG

Immunogen Information

Background: The protein encoded by this gene is a basic helix-loop-helix type transcription factor and shares similarity with Twist. This protein may inhibit osteoblast maturation and maintain cells in a preosteoblast phenotype during osteoblast development. This gene may be upregulated in certain cancers. Mutations in this gene cause focal facial dermal dysplasia 3, Setleis type. Two transcript variants encoding the same protein have been found.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: AMS; FFDD3; BBR SAY; Dermo1; SETLSS; bHLHa39; TWIST2

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human TWIST2 (NP_476527.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.