

NR5A2 Rabbit Polyclonal Antibody



CAB5766

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

43-70kDa

Calculated MW:

42kDa/53kDa/56kDa/61kDa

Applications:

WB IF

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IF 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

The protein encoded by this gene is a DNA-binding zinc finger transcription factor and is a member of the fushi tarazu factor-1 subfamily of orphan nuclear receptors. The encoded protein is involved in the expression of genes for hepatitis B virus and cholesterol biosynthesis, and may be an important regulator of embryonic development.

Immunogen information

Gene ID:

2494

Uniprot

O00482

Synonyms:

NR5A2; B1F; B1F2; CPF; FTF; FTZ-F1; FTZ-F1beta; LRH-1; LRH1; hB1F-2


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 206-495 of human NR5A2 (NP_003813.1).


Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot analysis of extracts of various cell lines, using NR5A2 antibody (CAB5766) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (CABS014) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.



Immunofluorescence analysis of A549 cells using NR5A2 antibody (CAB5766). Blue: DAPI for nuclear staining.