

Product Information

Product SKU:	CAB5797	Gene ID:	2118	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	54kDa	Conjugate:	Unconjugated
Calculated MW:	54kDa	Isotype:	IgG

Immunogen Information

Background:	Enables DNA-binding transcription activator activity, RNA polymerase II-specific and RNA polymerase II cis-regulatory region sequence-specific DNA binding activity. Involved in positive regulation of keratinocyte differentiation and positive regulation of transcription by RNA polymerase II. Located in chromosome and nucleolus.
Recommended Dilution:	WB,1:500 - 1:1000 IHC-P,1:50 - 1:200 IF/ICC,1:50 - 1:200
Synonyms:	E1AF; PEA3; E1A-F; PEAS3; Pea3 / ETV4
Purification Method:	Affinity purification
Immunogen:	Recombinant Protein corresponding to a sequence within amino acids 1-284 of human ETV4 (NP_001977.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.