

CAB5816

Product Information

Product SKU:	CAB5816	Gene ID:	10553	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	27kDa	Conjugate:	Unconjugated
Calculated MW:	27kDa	Isotype:	IgG

Immunogen Information

**Background:** Enables protein serine/threonine kinase activity. Involved in import into nucleus and regulation of angiogenesis. Acts upstream of or within positive regulation of programmed cell death; positive regulation of transcription by RNA polymerase II; and protein autophosphorylation. Located in cytosol and nuclear envelope.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	CC3; TIP30; SDR44U1; HTATIP2
Purification Method:	Affinity purification
Immunogen:	A synthetic peptide corresponding to a sequence within amino acids 100 to the C-terminus of human HTATIP2 (NP_006401.3).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.