

[KO Validated] HELLS Rabbit Polyclonal Antibody



CAB5831

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

120kDa

Calculated MW:

14kDa/36-39kDa/81-97kDa

Applications:

WB IHC IP

Reactivity:

Human, Mouse, Rat

Protein Background

This gene encodes a lymphoid-specific helicase. Other helicases function in processes involving DNA strand separation, including replication, repair, recombination, and transcription. This protein is thought to be involved with cellular proliferation and may play a role in leukemogenesis. Alternatively spliced transcript variants encoding different isoforms have been identified.

Immunogen information

Gene ID:

3070

Uniprot

Q9NRZ9

Synonyms:

HELLS; ICF4; LSH; Nbla10143; PASG; SMARCA6

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IHC 1:50
- 1:200 IP 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 719-838 of human HELLS (NP_060533.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images

