

CAB5842

---

## Product Information

<b>Product SKU:</b>	CAB5842	<b>Gene ID:</b>	4916	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

---

## Additional Information

<b>Observed MW:</b>	110-115kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	94kDa	<b>Isotype:</b>	IgG

---

## Immunogen Information

**Background:** This gene encodes a member of the neurotrophic tyrosine receptor kinase (NTRK) family. This kinase is a membrane-bound receptor that, upon neurotrophin binding, phosphorylates itself and members of the MAPK pathway. Signalling through this kinase leads to cell differentiation and may play a role in the development of proprioceptive neurons that sense body position. Mutations in this gene have been associated with medulloblastomas, secretory breast carcinomas and other cancers. Several transcript variants encoding different isoforms have been found for this gene.

**Recommended Dilution:** WB, 1:500 - 1:2000 IF/ICC, 1:50 - 1:200

**Synonyms:** TRKC; GP145-TrkC; gp145(trkC); TrKC

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 40-300 of human TrKC (NP\_001007157.1).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.