

Product Information

| | | | | | |
|--------------|---------|---------------|--------|-------------|-----------------|
| Product SKU: | CAB5895 | Gene ID: | 1106 | Size: | 20uL, 100uL |
| Clone No: | - | Host Species: | Rabbit | Reactivity: | Human,Mouse,Rat |

Additional Information

| | | | |
|----------------|--------|------------|--------------|
| Observed MW: | 260kDa | Conjugate: | Unconjugated |
| Calculated MW: | 211kDa | Isotype: | IgG |

Immunogen Information

| | |
|-----------------------|--|
| Background: | The CHD family of proteins is characterized by the presence of chromo (chromatin organization modifier) domains and SNF2-related helicase/ATPase domains. CHD genes alter gene expression possibly by modification of chromatin structure thus altering access of the transcriptional apparatus to its chromosomal DNA template. Alternatively spliced transcript variants encoding distinct isoforms have been found for this gene. |
| Recommended Dilution: | WB,1:500 - 1:2000 |
| Synonyms: | EEOC; DEE94; CHD2 |
| Purification Method: | Affinity purification |
| Immunogen: | Recombinant fusion protein containing a sequence corresponding to amino acids 1-223 of human CHD2 (NP_001036037.1). |
| Storage: | Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3. |