

CAB5922

Product Information

Product SKU: CAB5922	Gene ID: 30849	Size: 20uL, 100uL
Clone No: ARC2090	Host Species: Rabbit	Reactivity: Human,Mouse,Rat

Additional Information

Observed MW: 153kDa	Conjugate: Unconjugated
Calculated MW: 153kDa	Isotype: IgG

Immunogen Information

Background: Predicted to enable protein serine/threonine kinase activity. Involved in positive regulation of phosphatidylinositol 3-kinase activity; receptor catabolic process; and regulation of cytokinesis. Located in late endosome and microtubule cytoskeleton.

Recommended Dilution: WB,1:500 - 1:2000 IF/ICC,1:50 - 1:200

Synonyms: p150; VPS15; PIK3R4/VPS15

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 1259-1358 of human PIK3R4/VPS15 (Q99570).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.