

CAB5937

Product Information

Product SKU:	CAB5937	Gene ID:	64065	Size:	20uL, 100uL
Clone No:	ARC2094	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	21kDa	Conjugate:	Unconjugated
Calculated MW:	21kDa	Isotype:	IgG

Immunogen Information

Background: Involved in activation of cysteine-type endopeptidase activity. Predicted to be located in plasma membrane. Predicted to be active in cell-cell junction. Implicated in erythrokeratoderma variabilis and mutilating palmoplantar keratoderma with periorificial keratotic plaques.

Recommended Dilution:	WB,1:500 - 1:2000 IHC-P,1:50 - 1:200 IF/ICC,1:50 - 1:200
Synonyms:	THW; KCP1; EKVP7; OLMS2; PIGPC1; KRTCAP1; dJ496H19.1; PERP
Purification Method:	Affinity purification
Immunogen:	A synthetic peptide corresponding to a sequence within amino acids 94-193 of human PERP (Q96FX8).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.