

CAB5952

Product Information

Product SKU:	CAB5952	Gene ID:	94241	Size:	20uL, 100uL
Clone No:	ARC2102	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	27kDa	Conjugate:	Unconjugated
Calculated MW:	27kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable antioxidant activity. Involved in autophagic cell death; positive regulation of autophagy; and positive regulation of transcription, DNA-templated. Located in autophagosome; cytosol; and nucleus.

Recommended Dilution: WB,1:500 - 1:1000 IHC-P,1:50 - 1:200 IF/ICC,1:50 - 1:200

Synonyms: SIP; Teap; p53DINP1; TP53DINP1; TP53INP1A; TP53INP1B; p53 DINP1

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 100-200 of human p53 DINP1 (Q96A56).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.