

Product Information

Product SKU:	CAB5975	Gene ID:	27037	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	70kDa	Conjugate:	Unconjugated
Calculated MW:	69kDa	Isotype:	IgG

Immunogen Information

**Background:** The protein encoded by this gene is of unknown function. However, it is orthologous to the mouse Trmt2a gene and contains an RNA methyltransferase domain. Expression of this gene varies during the cell cycle, with aberrant expression being a possible biomarker in certain breast cancers. Several transcript variants encoding two different isoforms have been found for this gene.

Recommended Dilution:	WB,1:500 - 1:1000 IHC-P,1:50 - 1:200
Synonyms:	HTF9C; TRMT2A/HTF9C
Purification Method:	Affinity purification
Immunogen:	A synthetic peptide corresponding to a sequence within amino acids 500 to the C-terminus of human TRMT2A/HTF9C (NP_073564.3).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.