

Product Information

Product SKU:	CAB5977	Gene ID:	27336	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	129kDa	Conjugate:	Unconjugated
Calculated MW:	86kDa	Isotype:	IgG

Immunogen Information

**Background:** The protein encoded by this gene functions as a cofactor for the stimulation of transcriptional elongation by HIV-1 Tat, which binds to the HIV-1 promoter through Tat-TAR interaction. This protein may also serve as a dual-function factor to couple transcription and splicing and to facilitate their reciprocal activation. Alternatively spliced transcript variants have been found for this gene.

Recommended Dilution:	WB,1:1000 - 1:2000 IHC-P,1:50 - 1:100
Synonyms:	TATSF1; TAT-SF1; dJ196E23.2; HTATSF1
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human HTATSF1 (NP_055315.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.