

CAB6000

Product Information

Product SKU:	CAB6000	Gene ID:	81550	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	73kDa	Conjugate:	Unconjugated
Calculated MW:	73kDa	Isotype:	IgG

Immunogen Information

**Background:** Enables chromatin binding activity; methylated histone binding activity; and transcription coactivator activity. Predicted to be involved in chromatin organization and positive regulation of transcription, DNA-templated. Located in Golgi apparatus; cytosol; and nucleoplasm. Colocalizes with exon-exon junction complex.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	TDRD3
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human TDRD3 (NP_110421.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.