

UGT1A1 Rabbit Polyclonal Antibody

CAB6186



Product Information

Product SKU:	CAB6186	Gene ID:	-	Size:	20uL, 100uL
Clone No:	-	Host Species:	-	Reactivity:	-

Additional Information

Observed MW:	-	Conjugate:	-
Calculated MW:	-	Isotype:	-

Immunogen Information

Background:	-
-------------	---

Recommended Dilution:	-
Synonyms:	-
Purification Method:	-
Immunogen:	-
Storage:	-