

Product Information

Product SKU:	CAB6222	Gene ID:	57332	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	50kDa	Conjugate:	Unconjugated
Calculated MW:	43kDa	Isotype:	IgG

Immunogen Information

Background: Enables methylated histone binding activity. Involved in negative regulation of transcription by RNA polymerase II. Located in chromatin and nucleoplasm. Part of PRC1 complex. Biomarker of esophagus squamous cell carcinoma and glioblastoma.

Recommended Dilution:	WB,1:500 - 1:2000 IHC-P,1:50 - 1:200
Synonyms:	PC3; RC1; CBX8
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 80-280 of human CBX8 (NP_065700.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.