

ELF3 Rabbit Polyclonal Antibody



CAB6371

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

42kDa

Calculated MW:

39kDa/41kDa

Applications:

WB IHC IF

Reactivity:

Human, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IHC 1:50
- 1:200 IF 1:10 - 1:100

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

1999

Uniprot

P78545

Synonyms:

ELF3; EPR-1; ERT; ESE-1; ESX


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human ELF3 (NP_004424.3).


Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Product Images



Western blot analysis of extracts of various cell lines, using ELF3 antibody (CAB6371) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (CABS014) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (CABM00020). Exposure time: 90s.



Immunohistochemistry of paraffin-embedded rat kidney using ELF3 antibody (CAB6371) at dilution of 1:200 (40x lens).



Immunofluorescence analysis of A549 cells using ELF3 antibody (CAB6371). Blue: DAPI for nuclear staining.