

Product Information

Product SKU:	CAB6430	Gene ID:	8548	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	45kDa	Conjugate:	Unconjugated
Calculated MW:	45kDa	Isotype:	IgG

Immunogen Information

Background: Enables ubiquitin protein ligase binding activity. Acts upstream of or within Golgi organization and Golgi to plasma membrane protein transport. Located in Golgi apparatus and nucleoplasm. Biomarker of hepatocellular carcinoma.

Recommended Dilution:	WB,1:500 - 1:2000 IF/ICC,1:50 - 1:200
Synonyms:	JEM1; JEM-1; JEM-1s; GOLGIN-45; BLZF1
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-400 of human BLZF1 (NP_003657.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.