

NCR1 Rabbit Polyclonal Antibody



CAB6438

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to Figures

Calculated MW:

21kDa/22kDa/23kDa/32kDa/
34kDa

Applications:

WB

Reactivity:

Human

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

9437

Uniprot

O76036

Synonyms:

CD335; LY94; NK-p46; NKP46; NCR1

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 20-258 of human NCR1 (NP_001138929.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
