

PJA2 Rabbit Polyclonal Antibody



CAB6444

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to Figures

Calculated MW:

76kDa/78kDa

Applications:

WB IF

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IF 1:50 - 1:100

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

9867

Uniprot

O43164

Synonyms:

PJA2; Neurodap1; RNF131


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 439-708 of human PJA2 (NP_055634.3).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunofluorescence analysis of U2OS cells using PJA2 antibody (CAB6444). Blue: DAPI for nuclear staining.