

CAB6452

Product Information

Product SKU:	CAB6452	Gene ID:	10894	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	60kDa	Conjugate:	Unconjugated
Calculated MW:	35kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a type I integral membrane glycoprotein. The encoded protein acts as a receptor and binds to both soluble and immobilized hyaluronan. This protein may function in lymphatic hyaluronan transport and have a role in tumor metastasis.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: HAR; XLKD1; LYVE-1; CRSBP-1; LYVE1

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 20-238 of human LYVE1 (NP_006682.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.