

POMGNT2 Rabbit Polyclonal Antibody



CAB6484

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Calculated MW:

66kDa

Applications:

IF

Reactivity:

Human

Antibody Information

Recommended dilutions:

IF 1:20 - 1:50

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

This gene encodes a protein with glycosyltransferase activity although its function is not currently known.

Immunogen information

Gene ID:

84892

Uniprot

Q8NAT1

Synonyms:

POMGNT2; AGO61; C3orf39; GTDC2; MDDGA8


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 301-580 of human POMGNT2 (NP_116195.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunofluorescence analysis of A549 cells using POMGNT2 antibody (CAB6484). Blue: DAPI for nuclear staining.