

CAB6513

Product Information

Product SKU:	CAB6513	Gene ID:	2294	Size:	20uL, 100uL
Clone No:	ARC1863	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	40kDa	Conjugate:	Unconjugated
Calculated MW:	40kDa	Isotype:	IgG

Immunogen Information

Background: This gene belongs to the forkhead family of transcription factors which is characterized by a distinct forkhead domain. The specific function of this gene has not yet been determined; however, it may play a role in the regulation of pulmonary genes as well as embryonic development.

Recommended Dilution:	WB,1:500 - 1:1000
Synonyms:	FKHL5; ACDMPV; FREAC1; FOXF1
Purification Method:	Affinity purification
Immunogen:	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human FOXF1 (Q12946).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.