

CAB6756

Product Information

Product SKU:	CAB6756	Gene ID:	29966	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	87kDa	Conjugate:	Unconjugated
Calculated MW:	87kDa	Isotype:	IgG

Immunogen Information

Background: Enables armadillo repeat domain binding activity; protein phosphatase 2A binding activity; and small GTPase binding activity. Involved in negative regulation of transcription by RNA polymerase II and positive regulation of transcription by RNA polymerase II. Located in Golgi apparatus; nucleoplasm; and plasma membrane. Part of FAR/SIN/STRIPAK complex.

Recommended Dilution: WB,1:500 - 1:2000 IHC-P,1:100 - 1:200 IF/ICC,1:50 - 1:100

Synonyms: SG2NA; S/G2NA; PPP2R6B; STRN3

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 200-460 of human STRN3 (NP_001077362.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.