

# STRN3 Rabbit Polyclonal Antibody



CAB6756

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

87kDa

### Calculated MW:

77kDa/87kDa

### Applications:

WB IHC IF

### Reactivity:

Human, Mouse, Rat

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000 IHC

1:100 - 1:200 IF 1:50 -

1:100

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

29966

#### Uniprot

Q13033

#### Synonyms:

STRN3; PPP2R6B; S/G2NA; SG2NA

#### Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 200-460 of human STRN3 (NP\_001077362.1).

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---

