

Product Information

Product SKU:	CAB6772	Gene ID:	29765	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse,Rat

Additional Information

Observed MW:	40kDa	Conjugate:	Unconjugated
Calculated MW:	39kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable tropomyosin binding activity. Predicted to be involved in actin filament organization; muscle contraction; and myofibril assembly. Located in striated muscle thin filament.

Recommended Dilution:	WB,1:500 - 1:1000
Synonyms:	SK-TMOD; TMOD4
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-345 of human TMOD4 (NP_037485.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.