

TPPP3 Rabbit Polyclonal Antibody



CAB6775

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

19kDa

Calculated MW:

18kDa

Applications:

WB IHC

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IHC 1:50
- 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

51673

Uniprot

Q9BW30

Synonyms:

TPPP3; CGI-38; TPPP/p20; p20; p25gamma


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-176 of human TPPP3 (NP_057224.2).


Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Product Images



Western blot analysis of extracts of various cell lines, using TPPP3 antibody (CAB6775) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (CABS014) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (CABM00020). Exposure time: 300s.



Immunohistochemistry of paraffin-embedded human gastric cancer using TPPP3 antibody (CAB6775) at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded human colon carcinoma using TPPP3 antibody (CAB6775) at dilution of 1:100 (40x lens).