
Product Information

Product SKU:	CAB6793	Gene ID:	23390	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	73kDa	Conjugate:	Unconjugated
Calculated MW:	73kDa	Isotype:	IgG

Immunogen Information

Background: Enables identical protein binding activity and protein-cysteine S-palmitoyltransferase activity. Involved in lipoprotein transport and protein palmitoylation. Located in Golgi membrane; Golgi-associated vesicle membrane; and aggresome.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: HIP3; HYPH; HIP14; DHHC17; DHHC-17; HSPC294; ZDHHC17

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 170-310 of human ZDHHC17 (NP_056151.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.