

CSN1S1 Rabbit Polyclonal Antibody



CAB6892

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

20kDa/21kDa

Applications:

WB

Reactivity:

Mouse

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

1446

Uniprot

P47710

Synonyms:

CSN1S1; CASA; CSN1

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 16-135 of human CSN1S1 (NP_001881.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
