

CAB7114

Product Information

Product SKU:	CAB7114	Gene ID:	54957	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	17kDa	Conjugate:	Unconjugated
Calculated MW:	17kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to be involved in mRNA splicing, via spliceosome. Located in cytosol and nucleoplasm.

Recommended Dilution:	WB,1:500 - 1:2000 IF/ICC,1:50 - 1:100
Synonyms:	DLP; Dim2; TXNL4B
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-149 of human TXNL4B (NP_060323.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.